



# VENDYS® Digital Thermal Monitoring of Neurovascular Reactivity?



Recent unexpected findings showed that increasing fingertip temperature during cuff occlusion in the non-occluded (contra-lateral) arm may be associated with better vascular health, whereas decreasing temperature may be associated with higher cardiovascular risk. Upon confirmation, these findings may suggest a potential utility for **VENDYS®** to assess neurovascular as well as vascular reactivity.



Ischemia induced contra-lateral neurovascular vasodilation, mediated by nitric oxide?

